

**G85495**

## Spacesaver EHT Power Supply

**NFU 506**



### Purpose

For electrostatics experiments including ionisation, and for driving vacuum devices such as the Deflection e/m Tube.

### Output:

0 to 5000 volts dc, adjusted by rotary control, limited to 2mA, with a second output limited to 50 $\mu$ A, by a 50M $\Omega$  safety resistor. A 6.3V ac output is provided for the heater/filament of vacuum devices, maximum current 3A.

### Meter:

3 digit, backlit display of voltage, in kV.

### Ratings

230V ac +/-10% 50Hz 72VA

Fuse: T250mA time delay

Temperature range: 10 to 40 $^{\circ}$ C Humidity: 10 to 90% RH

Mass: 4.8kg

### EMC:

This equipment is Class A according to the EMC standard EN 55011 and is intended for use in a non-domestic environment only. Interference from nearby sources of RFI (radio frequency interference) may affect the accuracy of the display. Likewise, this power supply may interfere with nearby PCs, laptops, tablets and mobile phones. If problems are experienced, the items should be separated or switched off.

### Warnings

For your safety, this product should be used in accordance with these instructions, otherwise the protection provided may be impaired.

**Do not open or remove covers or panels.** Repairs and service may only be carried out by our repair agent, otherwise the warranty may be void.

Use only the 3-core mains cable supplied with the unit. If the mains cable is replaced, the rating of the replacement must be the same or better than the original.

The unit must be earthed at all times. The unit is earthed/grounded through the 3-core mains lead, so no additional earth connection is required.

Always position the power supply so that it can be disconnected from the mains, if an emergency arises.

This unit is intended for use in DRY conditions. Avoid spillage of water and other liquids on to the unit. If spillage occurs, disconnect the mains supply.

There is no specific requirement for insulation of external circuits as they cannot become hazardous live, as a result of connection to this unit. Limit the length of any connecting leads to 3 metres.

Shrouded plugs may be connected to the front panel sockets, to give extra protection from shocks, to users.

version G85495 2014



Caution, possibility of electric shock symbol

Earth (ground) terminal symbol

Caution

The high voltages produced by this unit are capable of giving shocks to careless users, but the maximum current is limited to 2mA, as recommended by CLEAPSS and SSERC. It is sensible for teachers and students with pacemakers or cochlear implants to take extra care, or avoid using the unit.

### Continuous use

This power supply is intended for school experiments in electrostatics, to power simple spectrum tubes, eg Geissler tubes, and for work with Teltron tubes where a low voltage heater/filament must also be powered. If the power supply is used continuously at full load, then we advise a maximum period of 60 minutes, followed by cooling for the same period. Similarly, if it is used at full load (3 amps) for 10 minutes, a cooling period of 10 minutes is advised.

### Features

- EHT output continuously variable by rotary control, 0 to 5kV dc at up to 2mA
- 6.3V ac output maximum current 3A
- 3 digit LCD indication of the output voltage
- Sockets accept standard 4mm plugs and shrouded plugs
- Illuminated mains on/off switch, prominently displayed.
- Internal fan for cool-running
- Time-delay fuse protection for transformer primary
- Additional thermal trip protection for transformer

### Outputs

The EHT and ac outputs may be used simultaneously, as required by the apparatus connected. The continuously variable 5kV output is fully isolated (floating).



The black socket used with the left-hand red socket, gives up to 2mA at 5kV.

The black socket used with the right-hand red socket, is limited to approximately 50 $\mu$ A at 5kV, by the internal safety resistor.

An earth terminal provides the facility to earth the negative or positive side of the output as required, allowing users to set up positive or negative potential differences with respect to earth (0V).

version G85495 2014

## Overload

If the EHT output is short circuited the voltage will “collapse” to zero and the short circuit should be removed as soon as possible. Certain vacuum tubes require more than 2mA to drive them. These cannot be used with this EHT supply. Attempts to draw more than 2mA result in the output voltage collapsing.

## Mains lead and fuse replacement


The mains cable stows on the back of the unit with a wind and clip feature. The mains cable retainer can be removed to allow the cable to be stored separately. Simply remove the screw and retainer from the base of the unit. The IEC cable is now detachable.



In the unlikely event of the fuse blowing, the retainer and mains cable must be removed to gain access to the fuse carrier.

The carrier can be removed with a fingernail or tip of a screwdriver.

A spare fuse is provided in the carrier.



Replacement fuses must be 20mm cartridge time-delay fuses, of the correct rating: T250mA

