

# Speed of sound

The Fast Timer does just what its name suggests!

It detects events, such as switches changing or light gates being interrupted, or in this case . . . loud sounds captured by microphones.

Times can be measured in milliseconds (0.001s) or microseconds (0.000001s).



## So let's do the speed of sound in air . . .

Connect the microphones to the sockets shown. Switch on the Fast Timer.

Check switch **a** is set to the right, and switch **b** is set to milliseconds.

Check that the microphones are 1 metre apart, and both are facing towards the hammer and steel plate.

Press **reset**.



Lift the steel plate off the table and strike it sharply with the hammer.

Check the time on the display. Press **reset** and repeat.

This is the time,  $t$ , in milliseconds, that sound takes to travel 1 metre.

Do the simple calculation 1000 divided by the time,  $t$

This is the speed of sound in metres per second.

## . . . and now in the table top

Ask someone to gently press the two microphones, face down, to the table top.

Check that they are 1 metre apart. Press **reset** on the Fast Timer.

Now hold the steel plate firmly on the table and strike it with the hammer.

You might need to set the Fast Timer to microseconds, to get a more accurate result.

The time is much shorter. Sound travels faster in solids (and liquids) than in air.

Can you calculate the speed of sound in the table top?

## Another possibility

If you had two plastic bags and an aquarium, you measure the speed of sound in water!