Important Safety Information

- Read the full manual before use. Safety may be impaired if it is not followed.
- Only connect to the voltage and frequency shown on the serial number plate. An earth is required.
- Always ensure the bath is filled with a liquid rated for the working temperature. See manual for liquid selection.
- Ensure the T100/TC120 is positioned so that the plug is easily accessible for disconnection.
- Before moving or cleaning, allow the unit to cool down and disconnect from the power supply.
- Safety can only be assured if the T100/TC120 is securely installed on a suitable vessel. See manual for guidance.


Features

1. Alarm Light
2. Display
3. F (function) Button
4. Pump Outlet Plate (*TC120 only*)
5. Heater On Light
6. Over-Temperature Dial (*TC120 only*)
7. S (select) Button
8. Dial
9. Float: low liquid level protection

Setting the temperature

1. Press S: Display will flash
2. Rotate dial to desired temp
3. Press S to store value


Setting the over-temp (TC120 only)

Rotate dial to desired temp cut-out (above set point)

Using the external pump (TC120 only)

Ensure the TC120 is switched off. Connect hoses securely to the pump outlet plate as shown left.

Always use pump connectors and hoses that are suitable for the operating temperature and liquid used.

Never disconnect any pipes or hoses while they contain very hot or very cold liquids or while the TC120 is pumping.

Safety cut out and alarm light

The T100/TC120 is fitted with a float switch which monitors the liquid level during operation. Should the liquid level drop below a certain level, the safety cut out will be activated, and alarm light will appear. Top up the bath with liquid until at least a minimum level to resume operation.

Cleaning

Clean the unit with a damp cloth, using only water. Do not use chemical cleaning agents.

Consult manual for full instructions.

Service & repair

For service or repair please contact our Service Department or your local distributor. All returned units must be accompanied by a Return Materials Authorisation (RMA) number, obtainable by contacting the Grant service department.

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