SFX/STEREO Series

## INSTRUCTION MANUAL

Model
MS-2
STX
SFX-31
SFX-32
SFX-33
SFX-34
SFX-51
SFX-52
SFX-91
SFX-91D
ST-30FX
ST-50Led

v 1.3 2017

## Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

## Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

### ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

## Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

## Intended use


For teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

---


## Overview MS-2


## Overview STX


## Overview SFX-31


## Overview SFX-32

## Overview SFX-33

## Overview SFX-34


## Overview SFX-51


## Overview SFX-52


## Overview SFX-91


## Overview SFX-91D

## Overview ST-30FX

## Overview ST-50Led

## Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

## Using the microscope

### 1. Place the specimen on the stage

Place the specimen on the microscope stage and lock the specimen using the slide-clamp if it is necessary. Ensure that the specimen is centred over the stage opening. (Fig.1)

Fig.1

### 2. Turn on the light

The microscope comes with an electrical illuminator. Insert the plug of the cable into the power socket, turn on the switch on the main body and select your light source. (Fig.2)

For SFX models, repeatedly press the touch button in order to change the light intensity. Depending on the model, the light intensity can be adjusted with a knob on the right side of the stand. (Fig.3-4)

Fig.2

Fig.3

Fig.4

### 3. Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. If two circles appear, the interpupillary distance is too big, and if two overlapped circles appear, the interpupillary distance is too small. (Fig.5)


Fig.5

#### 4. Focus

Rotate the focusing knob to bring the sample into focus with the highest magnification. (Fig.6)

Fig.6

#### 5. Dioptric compensation

This compensation makes it possible for people with glasses to adjust the microscope to their eyes and use the microscope without glasses. Adjust the diopter compensation ring of the right eyepiece tube until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece. (Fig.7)

Fig.7

#### 6. Magnification

Select the desired magnification by rotating the objective. (Fig.8)

Total magnification used can be calculated as:

Eyepiece magnification x Zoom magnification x Objective lens magnification.

Fig.8

#### 7. Contrast disc

You can use the black/white disc in order to enhance the contrast of the image when using the incident illumination.

If you want to look at bright opaque objects, place the disc with the black side up. (Fig.9)

Fig.9

---

## **Only for SFX-91D**

### **Installation of camera drivers and software**

Operating system requirements: Windows XP, Windows 7, Windows 8 or Windows 10. If your PC is Windows XP: before the installation of your Camera please check if Service Pack2 is installed. If not, please update your XP from the Microsoft web site. If your PC has a USB port less than 2.0, there may be some problems in capturing image from the camera due to the slow transmission.

1. Please insert your CD and run:

\optika vision lite \ setup.exe

2. Now Optika Vision Lite icon appears on your desktop.

3. Open the folder “\drivers” and run all the “.exe” files.

4. The installer may show warnings about missing Windows certification: ignore such warnings and select “Continue”. Note : do NOT cancel the installation.

5. At the end of the installation process restart the PC if requested.

6. Plug the Camera in your PC: the camera will be identified and Windows will recognize the driver which has just been installed from the CD.

Then ignore any message about Windows compatibility and select “Continue” .

Note : do NOT cancel the installation.

7. At the end of the identification process, disconnect your Camera and try to plug it again, in order to check if the Driver was installed correctly (no Windows error messages should appear).

Should any error message appear on your PC, you will find a “\Troubleshooting” folder on your CD which shows how to solve any problem which Windows could give.

### **Some notes**

You should know that your Windows PC needs to install a Driver before capturing images from the camera: after installing the driver from the CD, Windows needs to recognize the camera before using it.

Often a PC has some ports on the front and some on the back: the ports on the back are usually faster. We suggest to use these ones.

All our cameras, even though they may have an external connection, do not need any external power supply because it comes from the USB cable of your PC.

If you connect the camera into another USB port, Windows needs to recognize its Driver again even if it has already done this operation for the first USB port.

Generally, every time you use a peripheral device and after that you connect your Camera, then Windows needs to recognize the Driver again.

---

It may happen sometimes that Windows doesn't recognize a Camera using one of your PC USB ports, while gives no problem using another port.

To prevent any installation problem we have put on the CD all the software and drivers in specific folders with an appropriate names.

## **How to test your Camera**

In order to check if your Camera works correctly, it's not necessary to put it into the microscope: it's enough to connect it into the USB port, run Vision Lite using the icon on your desktop, select the right driver and check whether the shown image is not totally black.

## **Select the right drivers**

It's very easy:  
OPTIKA TCD 3.0

With Optika Vision Lite you have to tell the Software which is the Driver by selecting the microscope-shaped icon.

## **Getting started with Optika Vision Lite, a very easy measuring and documentation software.**

Run Optika Vision Lite by clicking on its icon on your desktop.

Select the right Driver for capturing the image. Note that after choosing the driver and before capturing the image you can set the brightness, saturation, gamma ecc..

Then you can capture the image.

The image is now ready for being manipulated using Optika Vision Lite.

# Maintenance

## Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

## To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

## Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

## Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.  
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## Equipment disposal

Art.13 DLsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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## **OPTIKA® S.r.l.**