



# ***Operating Instructions***

## ***5 Inch Digital Hotplate Stirrer***

***MI0102010***



# ***Contents***

<b>Safety Instructions</b>	<b>2</b>
<b>Inspection</b>	<b>4</b>
<b>Control</b>	<b>4</b>
<b>Trial run</b>	<b>6</b>
<b>Heating</b>	<b>6</b>
<b>Working with the external temperature sensor</b>	<b>7</b>
<b>Residual heat warning (Hot)</b>	<b>7</b>
<b>Stirring</b>	<b>8</b>
<b>Faults</b>	<b>8</b>
<b>Maintenance</b>	<b>8</b>
<b>Standards and regulations</b>	<b>9</b>
<b>Specifications</b>	<b>9</b>
<b>Warranty</b>	<b>10</b>

## Safety Instructions



For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions when in use.**
- Ensure that only trained staff work with the device.
- Socket must be earthed (protective ground contact)
- **Caution-Magnetism!** Effects of the magnetic field have to be taken in to account (e.g. data storage media, cardiac pacemakers...).
- The heating plate can reach temperatures of 310°C and pay attention to the residual heat after switching off.
- The main power supply cable should not touch the heating plate.
- Wear your personal protective equipment in accordance with the hazard category of the media to be processed. Otherwise there is a risk from:
  - Splashing and evaporation of liquids
  - Ejection of parts
  - Release of toxic or combustible gases
- Set up the device in a spacious area on an even, stable, clean, nonslip, dry and fireproof surface.
- Don't use damaged components.
- Gradually increase the speed and reduce the speed if
  - the medium splashes out of the container because the speed is too high
  - the machine is not running smoothly
  - the container moves on the heating plate
- Caution! Only process and heat up any media that has a flash point higher than the adjusted safe temperature limit that has been set.
- The safe temperature limit must always be set to at least 50°C lower than the fire point of the media used.
- Beware of hazards due to:

- flammable materials
- combustible media with a low boiling temperature
- glass breakage
- incorrect container use
- overfilling of media
- unsafe condition of container
- The machine may heat up when in use. Don't use the machine in explosive atmospheres with hazardous substances.
- Process pathogenic materials only in closed containers under a suitable extractor hood.
- Only process media that will not react dangerously to the extra energy produced in other ways, e.g. through light irradiation.
- The external temperature sensor PT100 must always be inserted in the media when connected and ensure it's inserted in the media to a depth of at least 20mm.
- Accessories must be securely attached to the machine and can't come off by themselves.
- Always disconnect the plug before fitting accessories.
- The machine can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300°C-400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.
- The voltage stated on the type plate must correspond to the mains voltage.
- Don't cover the machine, even partially e.g. with metallic plates or film, which will result in overheating.
- Ensure the heating plate is kept clean.
- Protect the machine and accessories from bumps and impacts.
- The minimum distance between the machines; the minimum

distance between the machine and the wall is mini. 100mm.

## Inspection

- Unpack the device carefully and check for any damages which may have arisen during transit. Please contact the manufacturer/supplier for technical support.

### Note:



If there is any apparent damage to the device, please do not connect the power line.



## Contents of package

Items	Qty
Main unit	1
Power cable	1
Screwdriver	1
Magnetic stirring bar	1
User Manual	1

## Control



Item	Description
① Knob for temperature setting	Set the temperature by rotating the knob
② Knob for speed setting	Set the speed by rotating the knob
③ Knob for safety temperature setting	Set the safety temperature by rotating the knob. When the heating temperature is higher than the safety temperature, the machine stops heating. It's the max. temperature of 320°C by rotating the knob to the limit in the clockwise direction; it'll be min. temperature of 50°C by rotating the knob to limit in anti clockwise direction.
④ Temperature setting scale	Reading the temperature setting value.
⑤ Speed setting scale	Read the speed setting value.
⑥ Heating indicator	It's lit when heating.
⑦ External temperature sensor indicator	When the external temperature sensor is connected, the indicator is lit and the LED displays the temperature of the external temperature sensor. When it's disconnected, the temperature of internal temperature sensor is displayed.
⑧ Setting indicator	The light will flash when setting the temperature or speed.
⑨ Power switch	Switch on/off

## ***Trial run***

- Make sure the required voltage corresponds to the mains voltage.
- Ensure the socket must be properly grounded.
- Switch on and start initializing.
- Add sample into the container with an appropriate stirring bar.
- Place the container on the working plate.
- Set the target speed and start stirring.
- Observe the stirring bar operation.
- Set the temperature and start heating
- Connect the external temperature sensor.

If the above operations are normal, the device is good to run. If these operations are abnormal, the device may probably be damaged during transit, please contact manufacturer/supplier for technical support.



### **Warning!**

DO NOT move the container when the device is working.

## ***Heating***

With the digital temperature control, the machine has two separate safe circuits. The hotplate is kept at a constant temperature by a digital control circuit. The hotplate temperature can also be monitored by another adjustable safe circuit. The two temperature sensors (PT1000)

- Set temperature by the temperature setting knob.
- Heating is switched on/off by pressing the temperature setting.
- The set temperature of last operation will be displayed when the machine is switched on. Normally, the set temperature and the actual temperature may probably have some differences:
  - Hotplate center and outer edge
  - The sample inside the container and the container.

To ensure the accuracy of the temperature inside the container,

please use the external temperature sensor PT100. internal for temperature control are built into the hotplate. The single external PT100 can monitor the temperature of sample.

### ***Working with the external temperature sensor***

The external temperature sensor PT100 is an optional accessory. When the sensor is connected, the sensor indicator will light and the machine is in the PT100 operation mode. The setting temperature of the external sensor and the actual temperature are displayed. The safety circuit controls the hotplate temperature.

Comparing with the temperature control of hotplate, the external temperature sensor can control the sample's temperature more precisely. The external temperature sensor must be inserted in the sample. For any faults, heating will be automatically stopped. In this case, please proceed as follows:

- Turn power off
- Ensure the external temperature sensor is inserted into the sample
- Turn on and set the temperature to start heating
- If the instrument still can't work properly, please contact manufacturer / supplier for technical support.

### ***Residual heat warning (HOT)***

To prevent the risk of burns from the hotplate, this instrument has residual heat warning function. When heating is stopped, and the heating plate temperature is still above 50°C, the setting indicator will light up to warn that there's potential risk of burns. When the hotplate temperature drops to below 50°C, the instrument will be automatically powered off. To turn off the instrument immediately, pull out the plug. In case of main power failure or disconnection, the residual heat warning will not work.



## **Stirring**

Stirring is switched on / off by pressing the speed setting knob. The speed is set by rotating the knob in the range of 50-1500rpm.

## **Faults**

- The device can't work when powered on
  - Check if the power line is connected properly
  - Check if the fuse is broken or loose
- Fault in power on self test
  - Power off the device and then restart
- Actual speed can't reach the set value
  - Media in high viscosity may probably cause abnormal speed reduction of the motor
- The device can't be powered off when switched off
  - Check if the residual heat warning function is still on and the hot plate temperature is above 50°C (the setting indicator is lit)

<b>Fault</b>	<b>Fault description</b>	<b>Solutions</b>
D1 light D3 flash	External temperature sensor PT100 failure	Replace PT100 and restart the device
D2 light D3 flash	Internal temperature sensor PT100 failure	Restart the device
D1/D2 light D3 flash	Motor failure	Reduce the sample

D1 for "Heating" indicator light, D2 for "External temperature sensor indicator light ; D3 for "Power on" indicator light.

If the instrument still can't work properly, please contact manufacturer / supplier for technical support.

## **Maintenance**

Proper maintenance can make the device work well and extend its life.

- Do not allow moisture to get into the device when cleaning.
- Disconnect the main plug when cleaning.
- Wear protective gloves when cleaning.
- Only use the recommended cleansing agents.

Dyes	isopropyl alcohol
Constructions materials	water containing surfactant / isopropyl alcohol
Cosmetics	water containing surfactant / isopropyl alcohol
Foodstuffs	water containing surfactant
Fuels	water containing surfactant

## ***Standards and regulations***

Construction in accordance with the following safety standards:

EN 61010-1

UL 3101-1

CAN/CSA C22.2(1010-1)

EN 61010-2-10

Construction in accordance with the following EMC standards:

EN 61326-1

Associated EU guidelines:

EMC-guidelines: 89/336/EWG

Instrument guidelines: 73/023/EWG

## ***Specifications***

Items	MI0102010
Voltage [VAC]	200-240/100-120
Frequency [Hz]	50/60
Power [W]	650
Stirring point position quantity	1
Max. stirring quantity (H <sub>2</sub> O) [L]	5
Max. magnetic bar [L x Φ, mm ]	55 x 10

Motor Type	Brushless DC motor
Max. power input of motor [W]	30
Max. power output of motor [W]	20
Speed range [rpm]	50-1500
Speed display	Scale reading
Temperature display	Scale reading
Working plate material	Aluminum alloy with ceramic coated
Dimension of working plate (mm)	Φ135
Heating power(W)	600
Temperature range[°C]	50-310
Safety temperature[°C]	50-310
External temperature sensor	PT100*
Control accuracy with the external temperature sensor [°C]	±1
Residual heat warning	50°C
Dimensions [W x D x H, mm]	250x148x85
Weight[kg]	1.9
Permissible ambient temperature[°C]	5-40
Permissible relative humidity	80%
Protection class acc. to DIN EN60529	IP42

\*Optional accessory

## **Warranty**

The instrument is warranted to be free from defects in materials and workmanship under normal use and service for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It doesn't cover any worn out parts, nor apply to any damage by improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

For claims under the warranty please contact your local supplier. You may also send the instrument directly to manufacturer, enclosing the invoice copy and giving reasons for the claim.