



# PH1001 Balance

[PH1001 – Revision A, Jan 2024)



## CONTENTS

1.0	INTRODUCTION .....	1
2.0	KEY OPERATIONS .....	1
3.0	CALIBRATION .....	2
4.0	MODEL SPECIFICATION .....	2
	WARRANTY INFORMATION .....	3

## 1.0 INTRODUCTION

Each PH1001 balance comes with 2 x AA batteries included. There is a choice of three capacities, each with three weighing units including grams, pounds and ounces, selectable by a Units key (grams only if the Mode key is not enabled). The bright backlit display enables operation in any degree of light, and the 304 grade stainless steel top pan offers protection against most materials and liquids. With easy calibration through the keypad the PH1001 ensures correct and accurate weighing at all times.

Always position the balance on a steady and level surface before use. Use the levelling feet and bubble indicator to adjust the balance to level.

The display “EEEE” indicates that the balance has been overloaded or damaged.

The display “Lo” or flashing battery indicator symbol indicates that the batteries must be changed.

Remove the cover on the underside of the balance to fit new batteries.

## 2.0 KEY OPERATIONS



Power on/off key.



Weighing mode selection key (if enabled) – g, oz or lbs.



Zero key is used to zero the balance or remove the weight of an empty container.

### 3.0 CALIBRATION

To calibrate the balance follow the instructions below.

- I. Turn on the scale and zero it.
- II. Press and hold **【MODE】** and **【→0←】** keys together for about 5 seconds until the display shows internal codes value.
- III. Press **【MODE】** key once, the display will show CALL→[xxxx], (where xxxx is the full capacity weight value).
- IV. Load calibration weight as per the displayed value, and the scale will be automatically calibrated. After calibration, the scale will return to normal weighing mode. Calibration is now complete.

### 4.0 MODEL SPECIFICATIONS

Model	PH1001
Maximum Capacity	1000g
Readability	0.1g
Units of Measure	grams (default), lbs and oz if the mode key is enabled
Operating Temperature	0 - 40 degrees C
Power Supply	2 x AA batteries or 6v DC 200ma Adaptor
Calibration Mass	1000g
Display	5 x 12mm/0.47 inches high digits
Pan Size	130mm/5.1 inches diameter
Overall Dimensions (w x d x h)	160 mm x 210mm x 65mm 6.3 in x 8.1 in x 2.6 in
Net Weight	385 g / 0.85 lbs

## WARRANTY INFORMATION

Limited Warranty (parts and labour) for the components failed due to defects in materials or workmanship. Warranty starts from the date of delivery.

During the warranty period, should any repairs be necessary, the purchaser must inform its supplier Phillip Harris. The company or its authorised Technician reserves the right to repair or replace the components at any of its workshops depending on the severity of the problems. However, any freight involved in sending the faulty units or parts to the service centre should be borne by the purchaser.


The warranty will cease to operate if the equipment is not returned in the original packaging and with correct documentation for a claim to be processed. All claims are at the sole discretion of Phillip Harris.

This warranty does not cover equipment where defects or poor performance is due to misuse, accidental damage, exposure to radioactive or corrosive materials, negligence, faulty installation, unauthorized modifications or attempted repair or failure to observe the requirements and recommendations as given in this User Manual. Additionally rechargeable batteries (where supplied) are not covered under warranty.

Repairs carried out under the warranty does not extend the warranty period. Components removed during the warranty repairs become the company property.

The statutory right of the purchaser is not affected by this warranty. The terms of this warranty is governed by the UK law. For complete details on Warranty Information, see the terms and conditions of sale available on our web-site.

## WEEE 2012/19/EU



This device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Disposal of batteries (if fitted) must conform to local laws and restrictions.  
Cet appareil ne peut être éliminé avec les déchets ménagers. L'élimination de la batterie doit être effectuée conformément aux lois et restrictions locales.  
Dieses Gerät nicht mit dem Hausmüll entsorgt.  
Dispositivo no puede ser desechado junto con los residuos domésticos  
Dispositivo non può essere smaltito nei rifiuti domestici.

### FCC / IC CLASS A DIGITAL DEVICE EMC VERIFICATION STATEMENT

NOTE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules and Canadian ICES-003/NMB-003 regulation. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

### CALIFORNIA PROPOSITION 65 - MANDATORY STATEMENT

WARNING: This product includes a sealed lead-acid battery which contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.



Products have been tested with, and are always supplied with mains power adaptors which meet all legal requirements for the intended country or region of operation, including electrical safety, interference and energy efficiency. As we often update adaptor products to meet changing legislation it is not possible to refer to the exact model in this manual. Please contact us if you need specifications or safety information for your particular item. Do not attempt to connect or use an adaptor not supplied by us.

ISO 9001:2008 certified global company with more than 40 years' experience in the production and sale of electronic weighing equipment.

All information contained within this publication is to the best of our knowledge timely, complete and accurate when issued. However, we are not responsible for misinterpretations which may result from the reading of this material.