

B8A23504/B8A23516/B8A23486

Murashige and Skoog

NFU 1259

Instructions

1. The required quantity should be mixed with about 75% of the total volume of water to be used, and then stirred to dissolve the powder.
2. Adjust the pH with 1M NaOH or HCl to 5.0+/-0.1 for liquid media. For media containing agar, adjust the pH to 5.7+/-0.1.
3. Then add the rest of the water to the final volume (plus agar if required). Heat and stir solution until clear. Other ingredients can be added at this stage, before dispensing into suitable containers and autoclaving at 15psi for 15 minutes.

Supplier details

Philip Harris Education, 2 Gregory Street, Hyde, Cheshire SK14 4RH

Orders and Information Tel: 0845 120 4521

Fax: 0800 138 8881

Repairs Tel: 0845 120 3211

E-mail: techsupport@philipharris.co.uk

Website: www.philipharris.co.uk

© Philip Harris Education, 2002, 2016

version 20.09