




Foundation Physics



Universal Stand

Universal Stand

KIT CONTENTS

ITEM	QUANTITY	DESCRIPTION
1	1	Universal Stand

OVERVIEW / SUGGESTIONS FOR USE

The UNIVERSAL Stand is designed to go with the Foundation Physics Line. It is designed to support each product and should not be used for any other products. It comes with adjustable screw legs and level so experiments can be performed on a level plane. The reservoir on the base is to hold extra marbles when not being used in an experiment.

ASSEMBLY

Insert the rod in the top hole and place the bolt up from the bottom as shown in the figure below. Turn the support stand clock wise until it is tight against the bottom of the base and the bolt is not sticking out form the bottom. The legs on the bottom can be unscrewed to adjust the level of the base.

