



150250

Tork Couch Roll Advanced

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in a healthcare environment. Rolls can be used in Tork Couch Roll Dispenser.

- Perforation for easy tear off and consumption control
- Strong and easy to roll, for optimal protection
- Better resistance and higher comfort due to 2-ply

[See more](#)



150250

Tork Couch Roll Advanced

Product data

System C1	Roll length 54.5 m
Roll width 50 cm	Roll diameter 12 cm
Number of sheets 165	Sheet length 33 cm
Core inner diameter 5.1 cm	Ply 2
Number of sheets 165	Embossing Yes
Color White	

Product certificates



The mark of
responsible forestry

150250

Tork Couch Roll Advanced

Shipping data

	Consumer units	Transport unit	Pallet
Height	500 mm	500 mm	1161 mm
Width	120 mm	360 mm	1000 mm
Length	120 mm	360 mm	1200 mm
EAN	7310791230351	7310791230344	7322540904192
Gross weight	1025.3 g	9.3 kg	167.54 kg
Net weight	980.1 g	8.8 kg	158.78 kg
Packing material	Plastic	Plastic	
Pieces	1	9	162