



ANTIBACTERIAL SURFACE CLEANSER

Findel code
1443304 (1L)
1424281 (500ml)

**KILLS 99.9% OF BACTERIA
AND COVID-19 VIRUS***

* Covid-19 Virus = SARS-CoV-2. Use disinfectants safely. Always read the label and product information before use.



Dettol Antibacterial Surface Cleanser is a surface spray that kills 99.9% of bacteria & viruses, including E. coli, salmonella, MRSA and the flu virus. Dettol Antibacterial Surface Cleanser is proven to remove: Allergens - Pollen Particles, Dust Mites, Pet Dander. This effective antibacterial spray is a non bleach, no taint and no odour solution that can be used on food preparation areas, highchairs and many more surfaces around the home.**



Kills 99.9% of bacteria and viruses**



Non bleach



No taint



No taint & no odour



Suitable for use on electronics***

Key Features

- Kills 99.9% of bacteria and viruses** including Covid-19 virus*
- Non bleach
- No taint & no odour
- Suitable for use on electronics***

Technical Specification

Country of Origin	GB
Form of packaging	BOS - Bottle with Trigger Spray
SKU Code	3003911
SKU Description	Dettol Antibacterial Surface Cleanser Trigger 750ml
Net Content/Size	750ml
EAN	5011417561935
Case Qty (units)	6
Cases PER PALLET	78

EN Test Method

Tested Organism(s)		Contact Time
EN 13697	E. coli, K. pneumoniae, E. hirae, S. aureus, P. aeruginosa	5 minutes
EN 13697	C. albicans	15 minutes
EN 14476	Bovine corona virus (BCV) Vacciniavirus Respiratory Syncytial Virus A (HRSV) Influenza A Virus (Avian Flu, H1N1)	30 seconds 5 minutes 5 minutes 30 seconds
EN 1650	C. albicans	20 minutes
EN 1276	P. aeruginosa, S. aureus, E. coli, E. hirae, S. aureus MRSA, S. enterica, L. monocytogenes, C. jejuni	5 minutes

Safety Information

Read label before use. Keep out of reach of children. If medical advice is needed, have product container or label at hand. Wash hands thoroughly after handling.

Not suitable for: Windows, mirror, fabrics, painted or varnished surfaces. Do not use with detergents or other disinfectants. Avoid contact with plastic baby bottles, which can crack when steam has sterilised. Do not freeze.

Instructions

1. Turn nozzle to ON position.
2. Simply spray directly on to surface.
3. Wipe over with CLEAN damp cloth & ALLOW TO DRY.
4. No need to rinse.
5. ***With electric equipment, spray on to cloth before wiping.



Ingredients

Per 100g of product contains 0.096g Benzalkonium Chloride. Contains < 5% Non-ionic surfactants, Disinfectant, Perfume.

* Covid-19 Virus = SARS-CoV-2. ** Enveloped viruses. Use disinfectants safely. Always read the label and product information before use.